

The screenshot displays the Windense software interface with two reports open. The left report is titled 'Compliance Report' and shows a table of lab performance metrics. The right report is titled 'Exception Report By Control' and shows a detailed list of control results for various labs.

Compliance Report Data:

QCS	TSR	Film	Paper	Printer
Central	Kemp	90.48% (19/21)	80.95% (17/21)	<None>
Central	Kilmer	92.86% (39/42)	88.10% (37/42)	<None>
Central	Stone	97.14% (34/35)	91.43% (32/35)	<None>
Central	COMPLIANCE %	93.88% (92/98)	87.76% (86/98)	<None>
Northern	Harris	91.43% (32/35)	91.43% (32/35)	<None>
Northern	Smith	85.71% (30/35)	91.43% (32/35)	<None>
Northern	Thomas	92.86% (26/28)	96.43% (27/28)	<None>
Northern	COMPLIANCE %	89.80% (88/98)	92.86% (91/98)	<None>
Western	Davis	78.57% (44/56)	78.57% (44/56)	<None>
Western	Blake	97.14% (34/35)	100.00% (35/35)	<None>
Western	Foster	73.81% (31/42)	78.57% (33/42)	<None>
Western	COMPLIANCE %	81.95% (109/133)	84.21% (112/133)	<None>

Exception Report Data:

Lab	1	2	3	4	5	6	7	8	9	10
0000001035	01035:Dehuxe	7257			8005555555				01/26/05	11:24
	Paper	01/11/05			Control					
	Paper	01/10/05			Control					
	Paper	01/09/05			< No Results >					
	Paper	01/08/05			Control					
	Paper	01/07/05			Control					
	57	: 71.43% Non-compliance								
	Film	01/13/05			Control					
	Film	01/12/05			Control					
	Film	01/11/05			Control					
	Film	01/10/05			Control					
	Film	01/09/05			Control					
	Film	01/08/05			Control					
	Film	01/07/05			Control					
	77	: 100.00% Non-compliance								
0000002857	02857:Dehuxe	1246			8005551212					
	Paper	01/13/05			< No Results >					
	Paper	01/12/05			< No Results >					
	27	: 28.57% Non-compliance								
	Film	01/13/05			< No Results >					
	Film	01/12/05			< No Results >					
	27	: 28.57% Non-compliance								
0000002888	02888:Dehuxe	2541			8005551212					
	Paper				100% Compliance					
	Film				100% Compliance					
0000002908	02908:Dehuxe	8095			8005555555					
	Paper	01/13/05			< No Results >					

The Compliance Report shows lab performance against the local rules for frequency of analysis. Results can be displayed by lab, by technical support representative or by regional manager (not shown here). This is a great tool for continuous improvement.